

New England Biolabs Product Specification

Product Name: Φ X174 DNA-Hae III Digest
Catalog #: N3026S/L
Concentration: 1,000 μ g/ml
Unit Definition: N/A
Shelf Life: 24 months
Storage Temp: -20°C
Storage Conditions: 10 mM Tris-HCl (pH 8.0), 1 mM EDTA
Specification Version: PS-N3026S/L v1.0
Effective Date: 24 Jun 2014

Assay Name/Specification (minimum release criteria)

A260/A280 Assay - The ratio of UV absorption of Φ X174 DNA-Hae III Digest at 260 and 280 nm is between 1.8 and 2.0.

DNA Concentration (A260) - The concentration of Φ X174 DNA-Hae III Digest is between 1000 and 1050 μ g/ml as determined by UV absorption at 260 nm.

Electrophoretic Pattern (Marker) - The banding pattern of Φ X174 DNA-Hae III Digest on a 3% agarose gel shows discrete, clearly identifiable bands at each band of the marker, when stained with Ethidium Bromide at a concentration of 0.5 μ g/ml.

Non-Specific DNase Activity (DNA, 16 hour) - A 50 μ l reaction in 1X NEBuffer 2 containing 5 μ g of Φ X174 DNA-Hae III Digest incubated for 16 hours at 37°C results in a DNA pattern free of detectable nuclease degradation as determined by agarose gel electrophoresis.



Derek Robinson
Director of Quality Control

Date 24 Jun 2014

