
New England Biolabs Product Specification

Product Name: BG-PEG-NH₂
Catalog #: S9150S
Shelf Life: 24 months
Storage Temp: -20°C
Specification Version: PS-S9150S v1.0
Effective Date: 09 Jun 2017

Assay Name/Specification (minimum release criteria)

Identity (Mass Spectrometry) - The observed molecular mass of BG-PEG-NH ₂ is 488.3 Da +/- 1 Da as determined by mass spectrometry analysis.

Physical Purity (HPLC) - BG-PEG-NH ₂ is ≥ 90% pure as determined by HPLC analysis.
--



Date 09 Jun 2017

Derek Robinson
Director of Quality Control

