

New England Biolabs Product Specification

Product Name: Monarch[®] Spin Columns S2D and Tubes
Catalog #: T1117L
Kit Components: Monarch[®] Spin Columns S2D (T1117)
Monarch[®] Spin Collection Tubes (T2118)

Shelf Life: 24 months
Storage Temp: 25°C
Specification Version: PS-T1117L v1.0
Effective Date: 25 Apr 2024

Assay Name/Specification (minimum release criteria)

Functional Testing (DNA Recovery and Purity) - Twenty-four Monarch [®] Spin Columns S2D are tested with 2 µg of Lambda DNA-BstEII Digest resulting in ≥75% recovery in ≥80% of the samples. OD 260/280 and 260/230 are ≥1.75 in ≥80% of the samples.
--

One or more products referenced in this document may be covered by a 3rd-party trademark.
Please visit www.neb.com/trademarks for additional information.



Date 25 Apr 2024

Lauren Brown
Quality Approver

